

Multiplex protein detection in 2-D gel electrophoresis using the Amersham ECL Plex fluorescent Western blotting system

Key words: *ECL Plex • fluorescent Western blotting • multiplex detection • two-dimensional gel electrophoresis • post-translational modification*

We examined the performance of the ECL Plex™ Western blotting system in both 1-D and 2-D Western blotting applications. The 1-D Western blotting experiment showed that both the nonphosphorylated and phosphorylated forms of GSK3β could be detected in PC-3U cells without TGF-β activation. The 2-D Western blotting experiment gave information on phosphorylation states that could not be obtained by 1-D Western blotting. A 48-kDa band corresponding to phosphorylated GSK3β in the 1-D experiments was resolved into at least five distinct protein isoforms, two of which were phosphorylated at serine 9.

Introduction

The Amersham™ ECL Plex Western blotting system based on sensitive CyDye™ conjugated secondary antibodies is most commonly used after separation of proteins on one-dimensional (1-D) SDS-PAGE gels. ECL Plex Western blotting provides sensitivity, linearity, and a dynamic range of nearly four orders of magnitude (1), but also the possibility of multiplex analysis—detection of two protein epitopes on the same membrane simultaneously. Multiplexing is commonly used for the quantitation of one protein relative to a protein of known abundance (housekeeping protein), but can also be used for the detection and quantitation of post-translational modifications such as phosphorylation, provided that there are antibodies available.

In this study, a two-dimensional (2-D) Western blotting approach was evaluated in addition to the traditional 1-D workflow, using fluorescent secondary antibodies, namely the ECL Plex Western Blotting Detection System. Human PC-3U cell lysate was separated by 1-D and 2-D SDS-PAGE, followed by blotting and detection of the beta form of glycogen synthase kinase-3 (GSK3β). ECL Plex CyDye-conjugated

secondary antibodies were used for simultaneous detection of the nonphosphorylated and phosphorylated forms of GSK3β. The results show that the ECL Plex Western blotting system enables multiplexing of GSK3β and Ser 9-phosphorylated GSK3β on 1-D gels, and that the 2-D Western approach can give higher resolution information on different phosphorylation states.

Materials

Products used

ECL Plex goat-α-mouse IgG-Cy3, 150 µg	PA43009
ECL Plex goat-α-rabbit IgG-Cy5, 150 µg	PA45011
CyDye DIGE Fluor, Cy™2 minimal dye	RPK0272
Hybond-LFP™, 20 x 20 cm, 10 sheets	RPN2020LFP
ECL Plex Fluorescent Rainbow™ Markers, full range, 500 µl	RPN851
ECL Advance™ Blocking Reagent	RPN418
2-D Clean-Up Kit	80-6484-51
2-D Quant Kit	80-6483-56
Immobiline™ DryStrip pH 7-11 NL, 7 cm	17-6003-68
DeStreak™ Rehydration Solution	17-6003-19
IPG Buffer pH 7-11 NL	17-6004-39
IPG Buffer pH 3-10 NL	17-6000-88
Ettan™ IPGphor™ II IEF System	80-6505-03
miniVE Vertical Electrophoresis System	80-6418-77
EPS 301 Power Supply	18-1130-01
TE 22 Mini Tank Transfer Unit	80-6204-26
Typhoon™ 9410 and ImageQuant™ TL	63-0055-80
PlusOne™ Bromophenol Blue	17-1329-01
PlusOne DTT	17-1318-01
PlusOne Glycerol	17-1325-01
PlusOne Glycine	17-1323-01
PlusOne SDS	17-1313-01
PlusOne Tris	17-1321-01



Other materials required

GSK3 beta antibody [G8] (abcam)	ab2602
Phospho-GSK-3beta (Ser9) Antibody (Cell Signaling)	9336S
10x PBS (Medicago)	12-9423-5
Novex™ 12% Tris-glycine gel (Invitrogen)	EC60055BOX
Novex 12% Tris-glycine gel, 2D well	EC6006BOX
Methanol (Merck)	K33730207
Tween™ 20 (Merck)	8.22184.1000
99.8% anhydrous dimethylformamide (Sigma-Aldrich)	705-6
Lysine (Sigma-Aldrich)	L-5626

Human prostate carcinoma (PC-3U) cells, a subclone of PC-3 [human prostate cancer cells (nr 3); ATCC number CRL-1435] modified for enhanced TGF- β receptor production (2)

Methods

1-D gel electrophoresis

Human PC-3U cells were lysed in a RIPA buffer (50 mM Tris-HCl, pH 8.0, 150 mM NaCl, 10% glycerol, 1% NP40, 0.5% deoxycholate) and loaded onto Novex 12% Tris-glycine gels in a series of four-fold dilutions from 30 μ g to 0.47 μ g of total protein. Gel electrophoresis was performed for 2.5 h at 100 V using the miniVE Vertical Electrophoresis System.

2-D gel electrophoresis

PC-3U cell lysate was further prepared for 2-D electrophoresis by using the 2-D Clean-Up Kit for desalting and cleanup. Protein concentration was measured with the 2-D Quant Kit. PC-3U cell lysate corresponding to 50 μ g total protein was labeled with 400 pmol of CyDye DIGE Fluor, Cy2 minimal dye for 30 min on ice. The reaction was stopped by the addition of 10 mM lysine.

For each Immobililine DryStrip gel, unlabeled PC-3U cell lysate (25 μ g of protein) was mixed with Cy2-labeled PC-3U cell lysate (5 μ g of protein).

First-dimension electrophoresis was performed on 7-cm Immobililine DryStrip pH 7-11 gels, selected because the predicted pI of GSK3 β is 8.95. A mixture of labeled and unlabeled human PC-3U cell lysate (30 μ g of protein) was applied to the rehydrated strips by cup loading, and the samples were focused for 11 kVh on the Ettan IPGphor II instrument (3).

Second-dimension SDS-PAGE was performed in the MiniVE Vertical Electrophoresis System using Novex 12% Tris-glycine gels with 2-D wells. ECL Plex Fluorescent Rainbow Markers were loaded onto the gel for size indication. Protein samples were first allowed to enter the gel gently for 15 min at 15 mA per gel, followed by separation for 2 h at 30 mA per gel.

Protein blotting and fluorescent detection

After electrophoretic separation the gels were transferred to Hybond-LFP membranes for 2.5 h at 25 V using a TE 22 Mini Tank Transfer Unit, followed by incubation in 2% ECL Advance Blocking Reagent in TBS overnight at 4°C. The blots were then incubated with rabbit polyclonal anti-phospho-GSK3 β (Ser9) and mouse monoclonal anti-GSK3 β primary antibodies (1:1000 dilution in TBS + 0.1% Tween-20; TBST) overnight at 4°C. The blots were washed twice quickly, then twice for 5 min each in TBST. They were then incubated for 1 h, protected from light, with the secondary antibodies: ECL Plex goat- α -rabbit IgG-Cy5 and ECL Plex goat- α -mouse IgG-Cy3 (1:2500 dilution in TBST). The membranes were washed three times quickly, then four times for 5 min each in TBST, followed by two brief washes in TBS. Membranes were dried at 40°C for 1 h before scanning.

Imaging was performed on the Typhoon 9410 scanner using the 633-nm (red) laser with a 670BP30 filter for Cy5-conjugated antibodies and the 532-nm (green) laser with a 580BP30 filter for Cy3 conjugates. The PMT value was adjusted until the most intense band nearly reached saturation. The images were then analyzed using ImageQuant TL software.

Application

Lysates of human PC-3U cells were prepared, run on 1-D and 2-D gels, transferred to Hybond-LFP membranes, and then probed with antibodies according to the protocol previously described. The targeted proteins were the non-phosphorylated and phosphorylated forms of GSK3 β (GSK3 β and pGSK3 β , respectively). As primary antibodies we used mouse monoclonal anti-GSK3 β targeted with ECL Plex goat- α -mouse IgG-Cy3 and rabbit polyclonal anti-phospho-GSK3 β (Ser9) targeted with ECL Plex goat- α -rabbit IgG-Cy5.

Results and discussion

One-dimensional gel electrophoresis

Using a mixture of ECL Plex goat- α -mouse IgG-Cy3 and ECL Plex goat- α -rabbit IgG-Cy5 fluorescent antibodies, GSK3 β and phospho-GSK3 β (Ser9) were simultaneously detected in one blot (Fig 1). The anti-GSK3 β primary antibody recognized two bands at about 48 and 51 kDa (Fig 1A). The 48-kDa band corresponds to GSK3 β and the 51-kDa band corresponds to GSK3 α , a result of cross-reactivity of the primary antibody, as stated by the manufacturer. The anti-pGSK3 β primary antibody also recognized two protein bands at 48 and 49-50 kDa (Fig 1B). The 48-kDa band corresponds to GSK3 β phosphorylated at serine 9. The 49- to 50-kDa band appeared somewhat fuzzy, not migrating as one distinct protein isoform. This may be because it corresponds to different phosphorylation states of GSK3 β , migrating differently in the gel. Another explanation could be that the 49- to 50-kDa band is nonspecific, that it does not correspond to GSK3 β .




Fig 1. 1-D electrophoresis of PC3U cell extract. A: ECL Plex Cy3-conjugated secondary antibody targeted against GSK3 β ; B: ECL Plex Cy5-conjugated secondary antibody targeted against the phosphorylated form of GSK3 β (pGSK3 β).

Figure 2 displays the same images as seen in Figure 1, but in color: images resulting from scanning with the green laser for GSK3 β (Fig 2A), the red laser for pGSK3 β (Fig 2B), and the overlay of scans from GSK3 β (Cy3, green) and pGSK3 β (Cy5, red) (Fig 2C). The yellow color of the 48-kDa band shows that the signals from the two primary antibodies overlap: the 48-kDa band contains phosphorylated GSK3 β .

TGF- β is a potent growth factor for stimulating a number of cellular responses, one of them resulting in the phosphorylation of GSK3 β (4). TGF- β activation will induce phosphorylation of GSK3 β significantly, but for this study only unactivated PC-3U cell lysate was used. We could not detect phosphorylation without any TGF- β activation; this can be explained by the ongoing dynamic process of phosphorylation/dephosphorylation in the cell. The overlay image of anti-GSK3 β and anti-phospho-GSK3 β (Ser9) scans clearly shows that a population of GSK3 β is phosphorylated in PC-3U cells without TGF- β activation.




Fig 2. Color images of 1-D Western blot scans from Fig 1. A: ECL Plex Cy3-conjugated secondary antibody targeted against GSK3 β ; B: ECL Plex Cy5-conjugated secondary antibody targeted against the phosphorylated form of GSK3 β ; C: overlay of A and B.

Two-dimensional gel electrophoresis

Two dimensional gel electrophoresis in mini-format allows high resolution, using both the pI and the molecular weight of proteins for separation, while the small format makes it convenient and fast compared to large-format 2-D electrophoresis. In this application, a fraction of the PC-3U cell lysate used for 2-D electrophoresis was prelabeled with CyDye DIGE Fluor Cy2 minimal dye to allow for direct visualization of total protein without post-staining of the membrane. Figure 3A shows the Cy2, Cy3, and Cy5 scans overlaid—total protein (blue), GSK3 β (green) and pGSK3 β (red). Figure 3B shows only the two forms of the target protein, GSK3 β (green) and pGSK3 β (red). Different epitopes at the same position detected by both antibodies appear yellow (Figs 3A and 3B). The positions of the detected spots agree with the predicted pIs of 8.95 and 8.98 and with the molecular weights of 48 and 51 kDa for the α and β forms of GSK3, respectively (http://au.expasy.org/tools/pi_tool.html).

Different phosphorylation states on 2-D gels are displayed as spot trains depending on the number of different states, where one phosphate will cause a shift towards the anode and a small increase in size (80 Da). GSK3 β has six potential phosphorylation sites (Ser9, Ser21, Thr43, Tyr216, Ser389, and Thr390; Cell Signaling). Using an isoelectric point calculator (5), the predicted pIs of GSK3 β nonphosphorylated and phosphorylated at one to six phosphorylation sites are 8.98 for nonphosphorylated GSK3 β , then 8.90, 8.81, 8.71, 8.59, 8.44, and finally 8.25 for the putative fully phosphorylated GSK3 β . It is likely that the spot trains seen in Figures 3A–C correspond to different phosphorylation states of GSK3 isoforms. The molecular weight of the two protein spots corresponding to phosphorylated GSK3 β in Figure 3 does not correspond to the molecular weight of the 49- to 50-kDa band seen in Figures 1B, 2B, and 2C, indicating that the 49- to 50-kDa band in the 1-D approach is not related to GSK3 β .

The results from the 2-D Western approach show that the 48-kDa band corresponding to phosphorylated GSK3 β in the 1-D experiments could be resolved into at least five distinct protein isoforms in the 2-D experiment, where two isoforms were phosphorylated at serine 9 (recognized by the anti-phospho-GSK3 β (Ser9) primary antibody).




Fig 3. Images of 2-D Western blot scans. A: Overlay of Cy2, Cy3, and Cy5 images; B: Overlay of Cy3 and Cy5 images; C: Cy3 image; D: Cy5 image. Total protein (Cy2, blue), ECL Plex Cy3-conjugated secondary antibody targeted against GSK3 β (green), and ECL Plex Cy5-conjugated secondary antibody targeted against the phosphorylated form, pGSK3 β (red). Molecular weights and pH range are shown in panel A. Protein spots representing GSK3 β (48 kDa), GSK3 α (51 kDa), and the serine 9-phosphorylated form of GSK3 β are indicated.

Conclusions

The high level of performance of the ECL Plex Western blotting system in both 1-D and 2-D Western blotting applications has clearly been demonstrated in this study. The 1-D Western blotting experiment showed that both the nonphosphorylated and phosphorylated forms of GSK3 β could be detected in PC-3U cells without TGF- β activation.

The 2-D Western blotting experiment gave additional information on phosphorylation states compared to 1-D Western blotting. The 48-kDa band corresponding to phosphorylated GSK3 β in the 1-D experiments could be resolved into at least five distinct protein isoforms, two of which were phosphorylated at serine 9.

Acknowledgment

The Isoelectric Point Calculator (PhosphoSite™, Cell Signaling) is based on an algorithm kindly provided by Elisabeth Gasteiger.

References

1. Application note: Multiplex protein detection using the ECL Plex fluorescent Western blotting system, GE Healthcare, 28-4015-40, Edition AA (2005).
2. Franzén, P. *et al.* Different signals mediate transforming growth factor-beta1-induced growth inhibition and extracellular matrix production in prostate carcinoma cells. *Exp. Cell Res.* **207**, 1–7 (1993).
3. 2-D Electrophoresis, Principles and Methods, GE Healthcare, 80-6429-60, Edition AC (2005).
4. Edlund, S. *et al.* Interaction between Smad7 and beta-catenin: importance for transforming growth factor beta-induced apoptosis. *Mol. Cell. Biol.* **25**, 1475–1488 (2005).
5. Bjellqvist, B. *et al.* The focusing positions of polypeptides in immobilized pH gradients can be predicted from their amino acid sequences. *Electrophoresis* **14**, 1023–1031 (1993).

Asia Pacific: Tel: +65 6275 1830 Fax: +65 6275 1829 **Australasia:** Tel: + 61 2 9899 0999 Fax: +61 2 9899 7511 **Austria:** Tel: 01/57606-1619 Fax: 01/57606-1627
Belgium: Tel: 0800 73 888 Fax: 02 416 82 06 **Canada:** Tel: 1 800 463 5800 Fax: 1 800 567 1008 **Central, East, & South East Europe:** Tel +43 1 972720
 Fax: +43 1 97272 2750 **Denmark:** Tel: 45 16 2400 Fax: 45 16 2424 **Eire:** Tel: 01494 544000 Fax: 0044 1494 542010 **Finland & Baltics:** Tel: +358-(0)9-512 39 40
 Fax: +358 (0)9 512 39 439 **France:** Tel: 01 6935 6700 Fax: 01 6941 9677 **Germany:** Tel: (089) 96281 660 Fax: (089) 96281 620 **Greater China:** Tel: +852 2100 6300
 Fax: +852 2100 6338 **Italy:** Tel: 02 27322 1 Fax: 02 27302 212 **Japan:** Tel: +81 3 5331 9336 Fax: +81 3 5331 9370 **Latin America:** Tel: +55 11 3933 7300
 Fax: +55 11 3933 7304 **Middle East & Africa:** Tel: +30 210 9600 687 Fax: +30 210 9600 693 **Netherlands:** Tel: 0800 82 82 82 1 Fax: 0800 82 82 82 4
Norway: Tel: 815 65 555 Fax: 815 65 666 **Portugal:** Tel: 21 417 7035 Fax: 21 417 3184 **Russia & other C.I.S. & N.I.S:** Tel: +7 (495) 956 5177 Fax: +7 (495) 956 5176
Spain: Tel: 93 594 49 50 Fax: 93 594 49 55 **Sweden:** Tel: 018 612 1990 Fax: 018 612 1910 **Switzerland:** Tel: 0848 8028 12 Fax: 0848 8028 13 **UK:** Tel: 0800 616928
 Fax: 0800 616927 **USA:** Tel: +1 800 526 3593 Fax: +1 877 295 8102

www.gehealthcare.com/ecl

GE Healthcare UK Limited
 Amersham Place
 Little Chalfont
 Buckinghamshire
 HP7 9NA
 UK



General Electric Company reserves the right, subject to any regulatory approval if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information. © 2006 General Electric Company—All rights reserved. GE and GE Monogram are trademarks of General Electric Company. Amersham, Cy, CyDye, DeStreak, ECL Advance, ECL Plex, Etan, Hybond, Hybond LFP, ImageQuant, Immobiline, PlusOne, Rainbow, and Typhoon are trademarks of GE Healthcare companies. Novex is a trademark of Invitrogen Corporation. PhosphoSite is a trademark of Cell Signaling Technology, Inc. Tween is a trademark of ICI Americas, Inc.

CyDye: this product or portions thereof is manufactured under license from Carnegie Mellon University under US patent number US5,268,486 and other patents pending. Some of these products may only be available to collaborators and customers within certain of our technology access programmes. The purchase of CyDye DIGE Fluors includes a limited license to use the CyDye Fluors for internal research and development, but not for any commercial purposes. A license to use the CyDye Fluors for commercial purposes is subject to a separate license agreement with GE Healthcare