

Mono Q PC 1.6/5

1 Introduction

The Precision Column PC 1.6/5 (0.10 ml) prepacked with Mono Q™ is designed for fast, sensitive and high resolution anion exchange separations of proteins, peptides and other biomolecules in the micropreparative scale.

The PC column is designed for use with the SMART™ system. To be able to use it with other high performance systems such as ÄKTAdesign™, you need the precision Column Holder (Code No 17-1455-01)

These instructions may be used as a guide for obtaining the best results from this column.

Additional information about column treatment, methods etc. can be found in the method manual for SMART system.

For further information about columns and equipment, please contact your local representative for GE Healthcare.

2 Important user information

- Filter buffers and samples to ensure long column life. (Section 5.21)
- Equilibrate the column before use. (Section 6.1)
- Use a flow rate within the recommended interval. (Section 6.2)
- Keep the column clean and store it properly (Sections 6.3, 7.1 and 7.2)

3 General description

3.1 Gel properties

Mono Q (Quarternary amino ethyl) is a strong anion exchanger based on a bed hydrophilic polymer. The monodispersity was accomplished through a unique process developed by Prof. John Ugelstad of SINTEF, Trondheim, Norway.

Mono Q has a particle size of 10 µm. The absence of fines

gives the packed column a large void volume (40%) and therefore a low back-pressure.

The charged group on the gel is -CH₂-N⁺-(CH₃)₃. Ionic capacity is 0.27–0.37 mmol/ml gel.

The maximum protein/peptide capacity is approximately 0.5–3 mg depending on the type of sample and running conditions. The practical loading range is 0.5–500 µg.

Nonspecific adsorption to MonoBeads™ is negligible. Recovery of protein activity is normally greater than 80%.

3.2 Chemical and physical stability

Mono Q PC 1.6/5 can be used with aqueous solutions in the pH range 2–12. It is resistant to solutions of urea up to 8 M, guanidine-HCl, ethylene glycol and similar compounds. The column is also stable in alcohol/water solutions.

Cationic or non-ionic detergents can be used with Mono Q. Do not use anionic buffers or detergents as they bind to the gel.

Dimethyl sulphoxide, dimethyl formamide, formic acid and similar solvents at concentrations <70% can be used with great care.

All oxidizing and other reactive substances should be avoided.

The maximum flow rate for Mono Q PC 1.6/5 is 0.4 ml/min, and the operational pressure limit is 5 MPa (50 bar, 750 psi). When running at a flow rate of 0.1 ml/min the column generally operates at a back pressure of 1.0–1.5 MPa.

All material in the column hardware is made of glass and fluoroplastics, whereas the filter is made of polypropylene and teflon. This composition makes the column totally biocompatible, and hence there is no risk of corrosion or sample contamination e.g. by metal ions.

The column operates at temperatures between +4 °C and +40 °C.

Table 1. Properties of Mono Q PC 1.6/5



Column dimensions	1.6 mm x 50 mm
Gel volume	0.10 ml
Charged group	-CH ₂ N ⁺ (CH ₃) ₃
Ionic capacity (CI)	0.27-0.37 mmol/ml gel
Average particle size	10 μm
Molecular weight range	10 ⁷
Maximum protein/peptide capacity	0.5-3 mg
Practical loading range	0.5-500 μg
Detection limit, one peak	10 ng/280 nm
Flow rate	0.01-0.4 ml/min
Sample volume	50 ml
Recovery of protein activity	80%
pH-stability	2-12
Operational pressure limit	5 MPa
Temperature interval	4-40 °C

4 Installation

4.1 Unpacking


Please check the contents against this list:

Description	No. per pack	Code No.
Mono Q PC I.6/5	1	17-0671-01
Instructions	1	
Quality control tests	1	

4.2 Connection of the column

The two protective plastic caps at either end of the column are held in place by bayonet fittings. To remove them follow the steps in figure 1.


Fig 1.



1. Hold the protective jacket tightly and press the cap inwards.
2. Turn the cap counter-clockwise.
3. Pull the cap outwards.


Follow the instructions in the software to install the column. Further information can be found in the system handbook. The procedure is illustrated in figures 2.1-2.4.

Fig 2.1. Check that the Automatic Column Holder is high enough to install



the column.

Fig 2.2. Put the column on to the Column Connector. The arrow on the label,



which should be pointing downwards, indicates the correct flow direction.

Fig 2.3. Align the column vertically with the Column Support





Fig 2.4. The column is now connected.



5. Getting started

5.1 Choice of running conditions

Anion exchangers should be used with cationic or zwitterionic buffers.

Cationic (e.g. cetyltrimethylammonium bromide) or non-ionic detergents (e.g. octylglucoside) may be used.

Avoid anionic buffers and detergents since they bind to Mono Q.

The pH of the buffer should be about 1 pH unit above the pI of the substances of interest.

Table 2. Suggested buffers for various pH intervals-

Operational pH intervals are approximately ± 0.5 pH units from pKa of the buffer. Based on experiments performed in our laboratories and determined at 25 °C.

Substance	pKa(25°C)	dpKa/°C ¹	Conc. (mM)	Counter-ion ²
N-methyl-piperazine	4.75	-0.015	20	Cl-
Piperazine	5.68	-0.015	20	Cl- /HCOObis-
Tris	6.46	-0.017	20	Clbis-
Tris propane	6.80		20	Cl-
Triethanolamine	7.76	-0.020	20	Cl-/CH ₃ COOTris
	8.06	-0.028	20	Cl-
N-methyl-diethanolamine	8.52	-0.028	50	Cl- /SO ₄ ²⁻ /CH ₃ COO
-1,3-diamino-propane	8.64	-0.031	20	Cl-
Ethanolamine	9.50	-0.029	20	Cl-
Glycine	9.90		20	Cl-
1.3-diamino-propane	10.47	-0.026	20	Cl-
Piperidine	11.12	-0.031	20	Cl-

1. The pKa value changes if the temperature is lowered.

2. Chloride is the most commonly used counter-ion, alternatives are listed in table 3.

Buffer concentration should be around 20–50 mM. For highly sensitive separations it is sometimes advisable to have a lower concentration, about 5–10 mM, to avoid buffer disturbances.

The UV-absorption profile of buffers can vary from batch to batch owing to contaminants. Perform a blank gradient run before running a sample at a high sensitivity.

The following is an example of a buffer system for anion exchange chromatography (provided the pI of the substance of interest is not too basic:

Start buffer: 20 mM Tris-HCl, pH 8–0

Elution buffer: 20 mM Tris-HCl, pH 8–0 + 0–5 M NaCl

Table 2 gives suggested buffers for various pH intervals.

Operational pH intervals are approximately ± 0.5

pH units from pKa of the buffer.

Specific anions have different elution strengths. Table 3 gives various anions and suggested concentrations.

Further information about chemical and physical stability is found in section 3.2.

Table 3. Anions and suggested concentrations for gradient elution.

Anion	Conc. (M)
O ²⁻	0.15
ClO ₄	0.15
I ⁻	0.20
BR ⁻	0.30
Cl ⁻	0.50
CH ₃ SO ₃ ⁻	0.60
HCOO ⁻	0.60
CH ₃ COO ⁻	0.70

Comment to table 3: Two millilitres is a reasonable gradient volume. Increase the anion concentration if the substance of interest is not eluted in the gradient.

5.2 Preparation of buffers and samples

The water for the buffers should be of ultrapure quality. Use HPLC grade solvents, salts and buffers. Degas and filter all solutions through a 0.22 µm filter. Be sure to select a solvent resistant filter if the buffer contains an organic solvent.

When possible, dissolve or dilute the sample in start buffer.

Pass the sample through a 0.22–0.45 µm filter. If necessary, centrifuge the sample first (10 000 × g for 10 min). Be sure to select a solvent resistant filter if the sample is dissolved in an organic solvent.

The sample may be applied directly to the column only if:

- the volume is very small and there is a risk of sample loss due to the dead-volume of the filter.
- the sample has previously passed through one or more chromatography steps, and is therefore estimated to be pure enough without filtration.

However, if a sample is applied directly there is a greater risk that the column will become blocked. Never apply a turbid solution (indicates sample insolubility.)

Note: Careful handling of solutions and samples increases column life considerably.

6. Running procedure

6.1 Column equilibration

The gel is delivered in 20% ethanol with chloride (Cl⁻) as the counter-ion. To equilibrate the column for first time use or after long term storage, proceed according to steps 1–4 below.

1. Install the column by following the instructions in section 4–2 and fill the system with the chosen buffers. Set the operational pressure limit to 5 MPa (50 bar, 750 psi).
2. Flush out the storage solution with 0.1 ml of start buffer at a flow rate of 0.1 ml/min.
3. Change to the desired counter-ion by washing with the eluent buffer. Continue until the system is completely free of air bubbles and the UV-signal is stable at the selected sensitivity. The minimum wash volume required is 0.5 ml.
4. Equilibrate with 0.5 ml of start buffer.

6.2 Sample application and elution

Ensure that the sample has been prepared as described in section 5.2.

The maximum protein/peptide capacity for Mono Q PC 1.6/5 is approximately 0.5–3 mg, but it is dependent on which proteins/peptides will be separated and which running conditions are used. Therefore a practical loading range is around 0.5–500 µg. This range minimizes the risk of losing non-bound material in the flow-through fractions, decreased resolution or detection problems.

The detection limit at 280 nm for a separate peak is about

10 ng. At 214 nm the detection limit is roughly 10 times lower.

The maximum sample volume is 50 ml. For volumes of 2 ml or less, use the sample loops listed in section 9.2. For volumes larger than 2 ml, use Superloop 10 ml or 50 ml (Code No. 19-7585-01 or 19-7850-01).

The recommended flow rate is 0.1 ml/min, but it can be varied between 0.01–0.4 ml/min.

Note: Depending on sample composition it may be necessary to reduce the flow rate during sample application to avoid a high back-pressure.

Use either linear or stepwise gradients for elution. A gradient volume of 2 ml (0.1 ml/min, 20 min) is usually sufficient. If the protein/peptide of interest is not eluted in the gradient, then increase gradient-time and/or anion concentration.

6.3 Column regeneration

To start a new run, clean the column with 0.5 ml of elution buffer, before equilibrating with 0.5 ml of start buffer. This procedure can be included in the method.

If it is suspected that not every thing has been eluted, inject 0.5 ml of a more concentrated salt solution (e.g. 1–2 M NaCl).

If detergents are not being used in the run, but have been included in previous ones, rinse the column with 0.5 ml of water before equilibrating with start buffer.

Should the column be dirty, then follow the cleaning procedures described in section 7.2 before equilibration is carried out.

7. Maintenance

7.1 Storage and prevention of microbial growth

Before storing for longer periods, wash the column sequentially with 0.5 ml of 1–2 M NaCl, 0.5 ml of water and 0.5 ml of 20% ethanol. Check that the UVsignal is stable before proceeding to the next wash.

To prevent draining during storage, always fill the empty space above the filters with liquid before replacing the protective caps. It is also advisable to keep the column wrapped in plastic or in a beaker filled with 200% ethanol.

Note: There is a risk that ethanol will leak into the space between the protective jacket and the column, but this will not affect the column.

The column should be stored between +4 °C and +30 °C. If stored for more than two weeks, place the column in a refrigerator.

7.2 Column cleaning

The following observations may indicate a dirty column which will require a more rigorous washing:

- increased back-pressure
- colour change in the gel bed
- loss of resolution
- decreased sample recoveries

Use own knowledge of possible contamination when choosing a cleaning method. The following is a list of recommended washing solutions. Inject them through a 2 ml loop (Code No. 18-5897-01) on to the column (flow rate <0.1 ml/min depending on the back pressure) and then rinse with 0.2–0.5 ml of water. Always ensure the baseline is stable before proceeding to the next wash.

- 2 ml 1 M NaCl
- 2 ml 2 M NaOH
- 2 ml 1 M HCl
- 2 ml 75% acetic acid
- 2 ml 1% TFA
- 1 ml 30% acetonitrile
- 2 ml 0.5% non-ionic detergent (e.g. octylglucoside) in 3 M HAC.

A common cleaning method used in our laboratories is:

- 1) 2 ml 2 M NaOH
- 2) 0.2 ml of water
- 3) 2 ml 1 M NaCl
- 4) 0.2 ml of water
- 5) 2 ml 75% HAC
- 6) 0.2 ml of water
- 7) 2 ml 1 M NaCl
- 8) 0.2 ml of water

If the column is still not restored, try injecting a solution of 1 mg/ml pepsin in 0.1 M acetic acid and 0.5 M NaCl and then leave it overnight at room temperature or 1 hour at 37 °C. Depending on the contamination other enzymes can be used, e.g. DNase. After the enzymatic treatment, repeat the cleaning process.

8. Trouble-shooting

8.1 Increased back-pressure

There are many reasons for increased back-pressure. One of the first checks should be the tubings and, if used, Prefilter 0.8 mm or 3.2 mm (Code No. 18-1800-75 and 18-1800-76).

Note: Prefilters are not meant to substitute the sample treatment described in section 5.2.

If increased back-pressure is caused by the column, try to restore it by performing the following steps :

1. Reverse the flow direction through the column by installing it upside-down and pumping 0.2 ml of elution buffer through it at 0.05 ml/min.
2. Return to normal flow direction and run for 5 min at 0.1 ml/min. If the problem persists, clean the column according to the description in section 7.2.

8.2 Checking the column packing

Though not as critical in gradient as in isocratic techniques, a well-packed column is important for a good separation. For best performance, ensure that the column is clean and has been correctly installed. If the column packing is the suspected cause of reduced resolution, repeat the function test described in figure 3 and compare the chromatograms.




Fig 3. Function test for Mono Q PC 1.6/5

9. Products and accessories

9.1 Precision Columns

The current range of Precision Columns available is listed below.

Column	Code No.
Ion exchange:	
Mini Q™ PC 3.2/3	17-0686-01
Mini S™ PC 3.2/3	17-0687-01
Mono Q PC 1.6/5	17-0671-01
Mono S™ PC 1.6/5	17-0672-01
Gel filtration:	
Superose™ 6 PC 3.2/30	17-0673-01
Superose 12 PC 3.2/30	17-0674-01
Superdex™ 75 PC 3.2/30	17-0771-01
Superdex 200 PC 3.2/30	17-1089-01
Superdex Peptide PC 3.2/30	17-1458-01
Reversed phase:	
µRPC C2/C18 SC 2.1/10	17-0704-01

9.2 Other products

Sample loops, teflon (5,50,100,200,500 µl)	18-0404-01
Sample loop 1000 µl + 2000 µl	18-5897-01
Union, M6 female/ 1/16" male, plastic	18-3858-01
Superloop 10 ml	19-7585-01
Superloop 50 ml	19-7850-01
Prefilter 0.8 mm	18-1800-75
Prefilter 3.2 mm	18-1800-76

9.3 Spare parts glass column

Protective caps (2/pk)	18-1015-28
------------------------	------------

www.gehealthcare.com/protein-purification
www.gehealthcare.com

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

Mono Q, SMART, ÄKTA, ÄKTAdesign, MonoBeads, Mini S, Mono S, Superose, Supersex and Drop Design are trademarks of GE Healthcare companies. GE, imagination at work and GE monogram are trademarks of General Electric Company

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

© 2006 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric company.

GE Healthcare Bio-Sciences AB
Björkgatan 30, SE-751 84 Uppsala, Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

Asia Pacific Tel: +85 65 62751830 Fax: +85 65 62751829 • Australasia Tel: +61 2 8820 8299 Fax: +61 2 8820 8200 • Austria Tel: 01 /57606 1613 Fax: 01 /57606 1614 • Belgium Tel: 0800 73 890 Fax: 02 416 8206 • Canada Tel: 1 800 463 5800 Fax: 1 800 567 1008 • Central, East, & South East Europe Tel: +43 1 972 720 Fax: +43 1 972 722 750 • Denmark Tel: +45 70 25 24 50 Fax: +45 45 16 2424 • Eire Tel: 1 800 709992 Fax: +44 1494 542010 • Finland & Baltics Tel: +358 9 512 3940 Fax: +358 9 512 39439 • France Tel: 01 69 35 67 00 Fax: 01 69 41 98 77 • Germany Tel: 0800 9080 711 Fax: 0800 9080 712 • Greater China Tel: +852 2100 6300 Fax: +852 2100 6338 • Italy Tel: 02 26001 320 Fax: 02 26001 399 • Japan Tel: 81 3 5331 9336 Fax: 81 3 5331 9370 • Korea Tel: 82 2 6201 3700 Fax: 82 2 6201 3803 • Latin America Tel: +55 11 3933 7300 Fax: +55 11 3933 7304 • Middle East & Africa Tel: +30 210 96 00 687 Fax: +30 210 96 00 693 • Netherlands Tel: 0800-82 82 82 1 Fax: 0800-82 82 82 4 • Norway Tel: +47 815 65 777 Fax: +47 815 65 666 • Portugal Tel: 21 417 7035 Fax: 21 417 3184 • Russia & other C.I.S. & N.I.S Tel: +7 495 956 5177 Fax: +7 495 956 5176 • Spain Tel: 902 11 72 65 Fax: 935 94 49 65 • Sweden Tel: 018 612 1900 Fax: 018 612 1910 • Switzerland Tel: 0848 8028 10 Fax: 0848 8028 11 • UK Tel: 0800 515 313 Fax: 0800 616 927 • USA Tel: +1 800 526 3593 Fax: +1 877 295 8102



imagination at work