

SOURCE

SOURCE 15PHE PE 4.6/100

Quick Information

SOURCE 15PHE PE 4.6/100 is a pre-packed PEEK column for rapid preparative hydrophobic interaction chromatography of proteins and oligonucleotides.

Read the Instruction

The instructions on this page will help you get started quickly with your new column. The other side gives more in-depth information on optimisation and trouble-shooting.

Column data

Matrix:	Polystyrene/divinyl benzene	
Ligand:	Phenyl	
Bead form:	Rigid, spherical porous, monodisperse	
Particle size:	15 µm	
Capacity:	At least 40 mg albumin/column	
pH stability:	Long term	Short term
	2–12	1–14
Temperature:	Regular use	Storage
	+4 to +40 °C	+4 to +30 °C
Pressure over the column	Regular use	Do not exceed
	0.25–2.5 MPa	4 MPa
Flow rate:	Regular use	Maximum
(water at room temperature)	0.5–2.5 ml/min	5.0 ml/min

First time use

Equilibration of the column before initial use or after long term storage or changing buffers:

1. 8 ml elution buffer.
2. 8 ml starting buffer.

Starting buffer: 50 mM phosphate buffer, 1.5 M (NH₄)₂SO₄, pH 7.0
 Elution buffer: 50 mM phosphate buffer, pH 7.0
 Flow rate: 1.0 ml/min

Note: Before connecting the column, start the pump and remove all air in the system, in particular in tubings and valves.



Try these conditions first

Before applying the sample, equilibrate the column. Proceed according to the section “First time use”.

Starting buffer*: 50 mM phosphate buffer, 1.5 M (NH₄)₂SO₄, pH 7.0
 Elution buffer: 50 mM phosphate buffer, pH 7.0
 Gradient: 0-100% B in 33 ml (20 CV)
 Flow rate: 2.0 ml/min

* Use a lower concentration of ammonium sulphate if the protein of interest begins to precipitate at this concentration.

Read the back of this instruction for information on optimising a separation.



Buffer and solvent resistance

De-gas and filter all buffers through a 0.22 µm filter.




Daily use

- Aqueous solutions pH 2–12
- Urea, up to 8 M
- Acetonitrile, up to 30% in aqueous buffers

Cleaning

- Acetonitrile, up to 30%
- Ethanol, up to 100%
- Methanol, up to 100%
- 2-propanol, up to 100%
- Hydrochloric acid, up to 1 M
- Sodium hydroxide, up to 2 M
- Acetic acid, up to 50%
- Guanidine hydrochloride, up to 8 M
- Anionic, cationic and non-ionic detergents



Avoid

- Unfiltered solutions
- Oxidising agents

Sample requirements/recommendations

Recommended ≤ 40 mg protein/column
 sample load:
 Sample preparation: Filtered through a 0.22 µm filter or centrifuged at 10 000 g for 10 min
 Temperature*: Ambient
 The sample should be dissolved in starting buffer.

*Hydrophobic interactions increases with increased temperature. Results achieved at room temperature may therefore not be reproduced in cold room, or vice versa.



In Depth Information

Delivery/storage

The gel is delivered in 20% ethanol. If the column is to be stored for more than two days after use, wash the column with at least 8 ml distilled water and then equilibrate it with at least 8 ml 20% ethanol.

Optimisation

Standard protocol

A first run can be performed as described in the section “Try these conditions first”. If the result is unsatisfactory, consider the following:

Action	Effect
Change pH/buffer*	Weaker/stronger binding
Change salt*	Selectivity change
Decrease the sample load	Improved resolution
Decrease the flow rate	Improved resolution
Use a shallower gradient	Improved resolution, but broader peaks and decreased concentration in fractions

* Refer to the section “choice of buffer”

Note: Hydrophobic interaction increases with increased temperature. Results achieved at room temperature may therefore not be reproduced in cold room, or vice versa.

Choice of buffers

In general, the adsorption process is often more selective than the desorption process and it is therefore important to optimize the starting “binding” buffer conditions with respect to:

- pH
- salt concentration
- type of salt

The combination of salt and pH can be manipulated to give optimum selectivity. Optimal conditions differ from application to application and are best established by running linear gradients and by varying the salt concentration and pH in the starting buffer.

The buffers given in the section “Try these conditions first” are recommended as the first choice.

The Hofmeister series, Table 1, gives approximate guidelines in choosing the type of salt to use. The most efficient salts are normally ammonium sulphate (up to 2 M) and sodium sulphate (up to 1 M) but also “weaker” salts such as sodium chloride (up to 4 M) can be considered.

Table 1. The Hofmeister series.

← Increasing binding effect

Anions: $\text{PO}_4^{3-} > \text{SO}_4^{2-} > \text{CH}_2\text{COO}^- > \text{Cl}^- > \text{Br}^- > \text{NO}_3^- > \text{ClO}_4^- > \text{I}^- > \text{SCN}^-$

Cations: $\text{NH}_4^+ > \text{Rb}^+ > \text{K}^+ > \text{Na}^+ > \text{Cs}^+ > \text{Li}^+ > \text{Mg}^{2+} > \text{Ca}^{2+} > \text{Ba}^{2+}$

← Increasing binding effect

For more information, please refer to the handbook “Hydrophobic Interaction Chromatography, Principles and Methods” available from Amersham Biosciences.

Column cleaning

Regular cleaning:

After each run, wash the column with 8 ml distilled water.

Before next run, re-equilibrate the column with at least 8 ml starting buffer until the UV base-line and pH/conductivity values are stable.

More rigorous cleaning:

Reverse the flow direction and run at a flow rate of 0.2 ml/min with the following sequence of solutions:

1. 5 ml 30% isopropanol
2. 5 ml distilled water
3. 5 ml 1 M NaOH
4. 5 ml distilled water.

Note: Do not store the column in 1 M NaOH.

Depending of the nature of the contaminants, the following cleaning solutions may also be appropriate:

- 70% ethanol
- 30% acetonitrile
- 2 M NaOH including 1 M NaCl
- 0.5% non-ionic detergent in 1 M acetic acid

Always rinse with 3 ml distilled water when any of the above cleaning solutions have been used. If detergents have been used, rinse with at

least 8 ml of 70% ethanol followed by 5 ml distilled water before equilibrating the column.

If column performance is still not restored, inject a solution of 1 mg pepsin/ml in 0.1 M acetic acid including 0.5 M NaCl and leave overnight at room temperature or 1 hour at 37 °C. Depending on the contamination, other enzymes can also be used, e.g. DNAase. After enzymic treatment, repeat step 1–4 in the rigorous cleaning procedure described above. Wash with elution buffer before equilibration with starting buffer and sample application.

Trouble-shooting

Symptom	Remedy
Increased back-pressure the column	Using reverse flow at 0.5 ml/min, pump 15 of ml elution buffer through the column. Then return to normal flow direction and run for 10 minutes at 2.0 ml/min.
Loss of resolution and/or decreased sample recovery	Follow the procedure described in the section “More rigorous cleaning”.
Air in the column	Reverse the flow direction and pump 20 ml of well de-gassed starting buffer at a flow rate of 0.5 ml/min.

DO NOT OPEN THE COLUMN!

Column performance control

We recommend checking the column performance at regular intervals. The function of the column can be checked as described in Figure 1.




Fig 1. Typical chromatogram from a function test of SOURCE 15PHE PE 4.6/100.

Accessories

Designation	No. per pack	Code No.
On-line Filter	1	18-1118-01
Handbook: “Hydrophobic Interaction Chromatography, Principles and Methods”	1	18-1020-90

Ordering information

Designation	No. per pack	Code No.
SOURCE 15PHE PE 4.6/100	1	17-5071-01

SOURCE and Drop Design are trademarks of Amersham Biosciences Limited or its subsidiaries.

Amersham is a trademark of Nycomed Amersham plc. Milli-Q is a trademark of Millipore Corp.

All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham Biosciences group which supplies them. A copy of these terms and conditions is available on request.

© Amersham Biosciences AB 2002 – all rights reserved

Amersham Biosciences AB
SE-751 84 Uppsala Sweden

Amersham Biosciences UK Limited
Amersham Place, Little Chalfont
Bucks HP7 9NA England

Amersham Biosciences Corporation
P.O. Box 1327
Piscataway, NJ 08855 USA

Amersham Biosciences Europe GmbH
Postfach 5480
D-79021 Freiburg Germany

Amersham Biosciences K.K.
Sanken Building, 3-25-1
Shinjuku-ku, Tokyo 169-0073 Japan