

Butyl Sepharose 4B

Butyl Sepharose™ 4B is a separation media for hydrophobic interaction chromatography (HIC). Substances are separated on the basis of their different hydrophobicities.

The butyl group is coupled to 4% agarose matrix by the cyanogen bromide method.



Table 1. Medium characteristics

Ligand density:	10–14 μmol butyl groups/ml drained medium
Available capacity:	Application dependent
Bead structure:	4% agarose
Bead size range:	45–165 μm
Mean particle size:	90 μm
Max linear flow rate*:	75 cm/h at 25 °C, HR 16/10 column, 5 cm bed height
pH stability**	
Long term:	4–8
Short term:	4–8
Chemical stability:	Stable to all commonly used aqueous buffers
Physical stability:	Negligible volume variation due to changes in pH or ionic strength
Autoclaving:	Not recommended

* Linear flow rate: $\frac{\text{volumetric flow rate (cm}^3\text{/h)}}{\text{column cross-sectional area (cm}^2\text{)}}$

** The ranges given are estimates based on our knowledge and experience. Please note the following:

pH stability, long term refers to the pH interval where the medium is stable over a long period of time without adverse effects on its subsequent chromatographic performance.

pH stability, short term refers to the pH interval for regeneration and cleaning.

Preparing the medium

Butyl Sepharose 4B is supplied pre-swollen in 20% ethanol. Prepare a slurry by decanting the ethanol solution and replace it with binding buffer, see below, in a ratio of 75% settled medium to 25% buffer. The binding buffer should not contain agents which significantly increase the viscosity. The column may be equilibrated with viscous buffers at reduced flow rates after packing is completed.

Packing Sepharose 4B

1. Equilibrate all material to the temperature at which the chromatography will be performed.
2. De-gas the medium slurry.
3. Eliminate air from the column dead spaces by flushing the end pieces with buffer. Make sure no air has been trapped under the column net. Close the column outlet with a few centimeters of buffer remaining in the column.
4. Pour the slurry into the column in one continuous motion. Pouring the slurry down a glass rod held against the wall of the column will minimize the introduction of air bubbles.
5. Immediately fill the remainder of the column with buffer, mount the column top piece onto the column and connect the column to a pump.
6. Open the bottom outlet of the column and set the pump to run at the desired flow rate. This should be at least 133% of the flow rate to be used during subsequent chromatographic procedures. However, the maximum flow rate, see Table 1, is typically employed during packing.
Note: If you have packed at the maximum linear flow rate, do not exceed 75% of this in subsequent chromatographic procedures.
7. Maintain the packing flow rate for 3 bed volumes after a constant bed height is reached.

Using an adaptor

Adaptors should be fitted as follows:

1. After the medium has been packed as described above, close the column outlet and remove the top piece from the column. Carefully fill the rest of the column with buffer to form an upward meniscus at the top.
2. Insert the adaptor at an angle into the column, ensuring that no air is trapped under the net.
3. Make all tubing connections at this stage. There must be a bubble-free liquid connection between the column and the pump, and column and the sample application system (LV-3 or LV-4 valves).
4. Slide the plunger slowly down the column so that the air above the net and in the capillary tubings is displaced by eluent. Valves on the inlet side of the column should be turned in all directions during this procedure to ensure that air is removed.
5. Lock the adaptor in position on the medium surface, open the column outlet and start the eluent flow. Pass eluent through the column at the packing flow rate until the medium bed is stable. Re-position the adaptor on the medium surface as necessary.

Equilibration

Before applying the sample, equilibrate the column with at least 2 column volumes at chosen binding conditions until the baseline is stable.

A common binding buffer for hydrophobic interaction chromatography is 50 mM phosphate buffer, 1.7 M $(\text{NH}_4)_2\text{SO}_4$, pH 7.0.

Binding

The binding of proteins to hydrophobic media is influenced by:

- The structure of the ligand *e.g.* carbon chain or an aromatic ligand. Butyl Sepharose 4B is less hydrophobic than Phenyl Sepharose CL-4B.
- The concentration and salting out effect of the binding buffer (see the Hofmeister series in Table 2).
Those salts which cause salting-out *e.g.* ammonium sulphate, also promote binding to hydrophobic ligands. The sample is applied in a solution of high concentration. A salt concentration between 0.5–2.0 M ammonium sulphate is commonly used*. The column should be equilibrated at the same concentration.
- Temperature. The hydrophobic interactions usually decrease with decreasing temperature.

Elution

Bound proteins are eluted by reducing the hydrophobic interaction. This can be done by:

- reducing the concentration of salting out ions in the buffer with a decreasing salt gradient (linear or step)
- eluting with a non-polar organic solvent, for example 40–50% ethylene glycol or 30% isopropanol
- including detergent in the eluent

* When working with proteins which have a tendency to aggregate, start with a lower $(\text{NH}_4)_2\text{SO}_4$ concentration to avoid the risk of precipitation.

Table 2. The Hofmeister series. Increasing the salting-out effect strengthens the hydrophobic interactions, whereas increasing the chaotropic effect weakens them.

	← Increasing salting-out effect								
Anions:	PO ₄ ³⁻	SO ₄ ²⁻	CH ₃ COO ⁻	Cl ⁻	Br ⁻	NO ₃ ⁻	ClO ₄ ⁻	I ⁻	SCN ⁻
Cations:	NH ₄ ⁺	Rb ⁺	K ⁺	Na ⁺	Cs ⁺	Li ⁺	Mg ²⁺	Ba ²⁺	
	Increasing chaotropic effect →								

As a start, we suggest a linear gradient from 0 to 100% B with:
 Buffer A: 50mM phosphate buffer, pH 7.0 + 1.7 M (NH₄)₂SO₄
 Buffer B: 50mM phosphate buffer, pH 7.0

Regeneration

Depending on the nature of the sample, regeneration is normally performed by washing with 2–3 bed volumes of an aqueous solution of 30% isopropanol and with 3 bed volumes of distilled water, followed by re-equilibrating in binding buffer.

Cleaning

In some applications, substances like denaturated proteins or lipids do not elute in the regeneration procedure. These can be removed by washing with 2 bed volumes of detergent in a basic or acidic solution. Use, for example, 0.1–0.5% non ionic detergent (*e.g.* Triton™ X-100) in 0.1 M acetic acid. Wash at a linear flow rate of approximately 40 cm/h, contact time 1–2 hours, reversed flow direction. When using high concentrations of organic solvents, apply increasing gradients to avoid air bubble formation.

Detergents are difficult to remove and are recommended only if the above procedures fail to give satisfactory results. After treatment with detergent wash with 5 bed volumes of 70% ethanol, or more, to remove residual detergents.

Alternatively, wash the column with 4 bed volumes of up to 70% ethanol.

Re-equilibrate the column with at least 3 bed volumes of binding buffer.

Storage

Butyl Sepharose 4B should be stored in neutral pH at 4–8 °C in presence of a bacteriostat, *e.g.* 20% ethanol.

Ordering Information

Designation	Pack size	Code No.
Butyl Sepharose 4B	50 ml	17-0960-01

Trademarks

Sepharose is a trademark of GE Healthcare Bio-Sciences Ltd.
GE Healthcare are trademarks of General Electric Company.
Drop Design is a trademark of GE Healthcare Bio-Sciences Ltd.
Triton is a trademark of Union Carbide
Chemicals and Plastic Company Inc.

All goods and services are sold subject to the terms and conditions of sale of the company within the General Electric Company group which supplies them. A copy of these terms and conditions is available on request.

GE Healthcare Bio-Sciences AB

Björkgatan 30
SE-751 84 Uppsala
Sweden

GE Healthcare UK Ltd

GE Healthcare Place, Little Chalfont
Buckinghamshire, England HP7 9NA

GE Healthcare Bio-Sciences Corp

800 Centennial Avenue
P.O. Box 1327
Piscataway, NJ 08855
USA

GE Healthcare Europe GmbH

Munzinger Strasse 9
D-79111 Freiburg
Germany

GE Healthcare Bio-Sciences KK

Sanken Building, 3-25-1
Shinjuku-ku, Tokyo 169-0073
Japan

