

GE Healthcare

Amersham
Megaprime™ DNA
Labelling Systems

Product Booklet

Codes: RPN1604
RPN1605
RPN1606
RPN1607



Page finder

1. Legal	3
2. Handling	4
2.1. Safety warnings and precautions	4
2.2. Storage and stability	4
2.3. Quality control	4
3. System components	6
3.1. Megaprime DNA labelling systems	8
4. Introduction	9
5. Megaprime DNA labelling protocols	11
5.1. Standard megaprime protocol	11
5.2. New megaprime protocol	15
5.3. Use of alternative reaction conditions	20
6. Appendices	27
6.1. Appendix I. Labelling of DNA fragments in low melting point agarose	27
6.2. Appendix II. Monitoring the reaction and calculating the specific activity of the labelled DNA	28
6.3. Appendix III. Removal of unincorporated nucleotides	32
6.4. Appendix IV. Additional equipment and reagents	34
7. Troubleshooting guide	35
8. References	38
9. Related products	39

1. Legal

GE and GE monogram are trademarks of General Electric Company. Amersham, Megaprime, Hybond, Hyperfilm, Hypercassette, Hyperscreen, Sensitize, Sephadex and SepRate are trademarks of GE Healthcare companies.

© 2006 General Electric Company – All rights reserved.

General Electric Company reserves the right, subject to any regulatory and contractual approval if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

Contact your GE Representative for the most current information and a copy of the terms and conditions

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited.
Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA UK

2. Handling

2.1. Safety warnings and precautions

Warning: For research use only.

Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

Caution: For use with radioactive material.

This product is to be used with radioactive material. Please follow the manufacturer's instructions relating to the handling, use, storage and disposal of such material.

All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls,

safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes, wash immediately with water. See material safety data sheet(s) and/or safety statement(s) for specific advice.

2.2. Storage and stability

Upon receipt of these systems components should be stored at -15°C to -30°C. The components are stable for at least 3 months when stored under recommended conditions.

2.3. Quality control

The Megaprime DNA labelling systems are tested by our quality control group to ensure an incorporation rate greater than 55% after 10 minutes at 37°C.

The performance of RPN 1604/1605 is tested with the standard DNA provided

using 17 pmol/25 ng DNA of [α - 32 P] labelled nucleotides, specific activity 3000 Ci/mmol (codes PB 10204-7) and RPN 1606/1607 are tested using 17 pmol/25 ng DNA of [α - 32 P]dCTP, 3000 Ci/mmol (code PB 10205). Incorporations greater than 55% are achieved after 10 minutes incubation at 37°C, as assayed by thin-layer chromatography on PEI cellulose in 1.25 M KH_2PO_4 . PH3.4.

In addition components of the kits are checked for identity by HPLC and the DNA solutions for concentration by UV spectrophotometry.

3. System components

Magaprime DNA labelling	RPN1604	RPN1605	RPN1606	RPN1607
Primer solution: Random nonamer primers in an aqueous solution	150 µl	300 µl	150 µl	300 µl
Labelling buffer; dATP, dGTP and dTTP in Tris/HCl pH7.5, 2-mercaptoethanol and MgCl ₂	-	-	300 µl	600 µl
Nucleotide solutions				
(a) dATP	120 µl	240 µl	-	-
(b) cCTP	120 µl	240 µl	-	-
(c) dGTP	120 µl	240 µl	-	-
(d) dTTP	120 µl	240 µl	-	-
in Tris/HCl pH8.0, 0.5 mM EDTA				
Reaction buffer: A 10x concentrated buffer containing Tris/HCl pH7.5, 2-mercaptoethanol and MgCl ₂	150 µl	300 µl	-	-

Magaprime DNA labelling	RPN1604	RPN1605	RPN1606	RPN1607
--------------------------------	----------------	----------------	----------------	----------------

Enzyme solution; 1 unit/ μ l DNA polymerase 1 Klenow fragment (cloned in 100 mM potassium phosphate pH6.5, 10 mM 2-mercapto- ethanol and 50% glycerol	60 μ l	120 μ l	60 μ l	120 μ l
---	------------	-------------	------------	-------------

Standard DNA solution; 5 ng/ μ l <i>Hind</i> III digested lambda DNA in 10 mM Tris/HCl pH 8.0, 1 mM EDTA	25 μ l	50 μ l	25 μ l	50 μ l
--	------------	------------	------------	------------

Carrier DNA solution; 500ng/ml sonicated herring sperm DNA in 10 mM Tris/HCl pH 8.0, 1 mM EDTA	1.25 ml	2.5 ml	1.25 ml	2.5 ml
--	---------	--------	---------	--------

3.1. Megaprime DNA labelling systems

30 standard labelling reactions –
for use with any radioactive nucleotide RPN 1604

60 standard labelling reactions –
for use with any radioactive nucleotide RPN 1605

30 standard labelling reactions –
for use with radioactively labelled dCTP RPN 1606

60 standard labelling reactions –
for use with radioactively labelled dCTP RPN 1607

4. Introduction

Feinberg and Vogelstein (1,2) introduced the use of random sequence hexanucleotides to prime DNA synthesis on denatured template DNA at numerous sites along its length. The primer-template complex is a substrate for the 'Klenow' fragment of DNA polymerase 1. By substituting a radiolabelled nucleotide for a non-radioactive equivalent in the reaction mixture newly synthesized DNA is made radioactive (see Figure 1). The absence of the 5'-3' exonuclease activity associated with DNA polymerase 1 ensures that labelled nucleotides incorporated by the polymerase are not subsequently removed as monophosphates. Very small amount of input DNA can be labelled, enabling very high specific activity DNA probes to be produced with relatively small quantities of added nucleotides. These radioactive labelled fragments can then be used as sensitive hybridization probes for a wide range of filter based applications (3-6).

Previous protocols for the random primer labelling of DNA have required reaction times of at least 30 minutes. GE Healthcare's Megaprime DNA labelling system allows the labelling of template DNA to the same high specific activity but at a greatly accelerated rate. Probes of specific activity 1.9×10^9 dpm/ μ g can be produced with the majority of DNA substrates, using the standard protocol, after 10 minutes incubation at 37°C. This rapid labelling is achieved by the use of nonamer primers rather than the conventional hexamers (Figure 1). Nonamers allow for more efficient priming from the template DNA at 37°C, resulting in fast and efficient labelling of the DNA. A new alternative protocol has further reduced the variability in labelling which can occur with DNA template from a variety of sources. Both the standard Megaprime protocol and the new protocol are given as options in this booklet. The labelling of DNA in low melting point agarose takes only 15-30 minutes in contrast to conventional systems where overnight incubation are necessary.




Figure 1. Preparation of labelled probes using GE Healthcare's megaprime DNA labelling systems.

5. Megaprime DNA labelling protocols

The Megaprime systems allow DNA from a variety of sources to be labelled *in vitro* to high specific activity with ^{32}P and other radionuclides. The specific activity of the probes generated by these systems will vary according to the specific activity of the labelled dNTP used.

The standard Megaprime protocol is presented, together with a new protocol which reduces the variation in labelling efficiency that can occur with DNA template from a variety of sources.

The protocols given here are for use with 17 pmol[α - ^{32}P]dNTP, specific activity 3000 Ci/mmol. For alternative reaction conditions refer to page 20.

DNA prepared by standard minilysate methods may be used in either protocol. DNA solutions which are too dilute to be used directly should be concentrated by ethanol precipitation followed by redissolution in an appropriate volume of water or 10 mM Tris/HCl, pH 8.0, 1 mM EDTA. DNA in restriction enzyme buffers may be added directly to the reaction. The reaction can also be performed with DNA in agarose gel slices (see note 3 and Appendix 1).

5.1. Standard Megaprime protocol

Protocol	Notes
1. Dissolve the DNA to be labelled to a concentration of 2.5–25 ng/ μl in either distilled water or 10 mM Tris/HCl, pH8.0, 1 mM EDTA (TE buffer).	1. If desired, the labelling efficiency of a DNA sample can be compared with that of the standard DNA supplied with the kit. In this case 5 μl of standard DNA should be used.

Protocol**Notes**

- Place the required tubes from the Megaprime system, with the exception of the enzyme, at room temperature to thaw. Leave the enzyme at -15°C to -30°C until required, and return immediately after use.
 - Place 25 ng of template DNA into a microcentrifuge tube and to it add 5 μl of primers and the appropriate volume of water to give a total volume of 50 μl in the final Megaprime reaction. Denature by heating to $95\text{--}100^{\circ}\text{C}$ for 5 minutes in a boiling water bath.
 - Spin briefly in a microcentrifuge to bring the contents to the bottom of the tube.
 - Keeping the tube at room temperature, add the nucleotides and reaction buffer (RPN 1604/5) or the labelling buffer (RPN 1606/7) followed by the radiolabelled dNTP(s) and enzyme as follows:
- When labelling DNA in low melting point agarose, first place the tube containing the stock DNA in a boiling water bath for 30 seconds to melt the agarose before removing the required volume. The volume of low melting point agarose DNA should not exceed 25 μl in a 50 μl reaction.
 - The reaction volume may be scaled up or down if more or less than 25 ng of DNA is to be labelled.

Protocol**Notes**

Component	RPN1604/5	RPN1606/7
Labelling buffer		10 μ l
Unlabelled dNTPs	4 μ l of each omitting those to be used as label	-
Reaction buffer	5 μ l	-
Radiolabelled (dNTP)	5 μ l	5 μ l (dCTP)
Enzyme	2 μ l	2 μ l

6. Mix gently by pipetting up and down and cap the tube. Spin for a few seconds in a microcentrifuge to bring the contents to the bottom of the tube.

7. Incubate at 37°C for 10 minutes

6. Avoid vigorous mixing of the reaction mixture as this can cause severe loss of enzyme activity.

7. Purified DNA can be labelled to high specific activity in 10 minutes at 37°C but, if desired, can be labelled for up to 1 hour at this temperature. When labelling DNA in low melting point agarose, longer incubation of 15–30 minutes at 37°C are required for optimum labelling. Longer incubation

Protocol

7. Incubate at 37°C for 10 minutes *continued*.
8. Stop the reaction by the addition of 5 µl of 0.2 M EDTA. For use in a hybridization, denature the labelled DNA by heating to 95–100°C for 5 minutes, then chill on ice.

Notes

7. **Continued.** times (up to 60 minutes) are required when nucleotide analogues (e.g. [³⁵S]dNTPαS) are used.
8. Labelled probe can be stored at -15°C to -30°C in a non frost-free freezer. Prolonged storage of ³²P-labelled probes can lead to substantial probe degradation(7). High specific activity probes should be stored for no longer than 3 days. Although probe purification is not usually necessary for most membrane applications, the removal of unincorporated nucleotide is sometimes useful to reduce background in filter hybridizations for probes >10⁹ dpm/µg or when the reaction yields an incorporation of less than 50%. This procedure is described in Appendix III. Calculation of probe specific activity is described in Appendix II. Extensive experimentation with **Rapid-hyb buffer (RPN1635/6)** has shown that probe purification, even

Protocol	Notes
<p>8. Stop the reaction by the addition of 5 μl of 0.2 M EDTA. For use in a hybridization, denature the labelled DNA by heating to 95–100°C for 5 minutes, then chill on ice <i>continued</i>.</p>	<p>8. Continued under the conditions given above is not required with the isotopes ^{32}P and ^{33}P. Purification of ^{35}S labelled probes is however required to reduce filter background.</p>

5.2. New Megaprime protocol

Protocol	Notes
<p>1. Dilute the DNA to a concentration of 5 ng/μl in either distilled water or 10 mM TE buffer.</p> <p>2. Place the required tubes from the Megaprime system with the exception of the enzyme at room temperature to thaw. Leave the enzyme at -15°C to -30°C until required, and return immediately after use.</p>	<p>1. DNA solutions at concentrations in the range 5–25 ng/μl can be used if desired. However the denaturing volume (step 3) should not be less than 10 μl to maximize the efficiency of primer annealing. The labelling efficiency of a DNA sample can be compared with that of the standard DNA supplied with the kit. In this case 5 μl of standard DNA should be used.</p>

Protocol

Notes

3. Place 25 ng (5 μ l) of template DNA into a clean microcentrifuge tube and to it add 5 μ l of primers. Denature by heating to 95–100°C for 5 minutes in a boiling water bath.
4. Spin briefly in a microcentrifuge to bring the contents to the bottom of the tube.
5. Keeping the tube at room temperature add the nucleotides and 10x reaction buffer (RPN 1604/5) or the labelling buffer (RPN 1606/7), water and enzyme:-

3. If the volume of DNA and primers is less than 10 μ l make up to this volume with water. When labelling DNA in low melting point agarose first place the tube containing the stock DNA in a boiling water bath for 30 seconds to melt the agarose before removing the required volume. The volume of low melting point agarose DNA should not exceed 25 μ l in a 50 μ l reaction.
5. The enzyme can be added directly to the reaction mix or pipetted on to the side of the microcentrifuge tube and “washed” down with the water.

Component	RPN1604/5	RPN1606/7
Labelling buffer		10 μ l
Unlabelled dNTPs	4 μ l of each – omitting those to be used as label	

Protocol**Notes**

Reaction buffer	5 μ l	-
Enzyme	2 μ l	2 μ l
Water*	as appropriate for a final reaction volume of 50 μ l*	

* When calculating this volume remember to allow for the volume of radioactive nucleotide to be added.

6. Cap the tube and spin for a few seconds in a microcentrifuge to bring the contents to the bottom of the tube.

7. Add the radiolabelled dNTP, for example 5 μ l [α -³²P]dNTP, specific activity 3000 Ci/mmol. Mix by gently pipetting up and down. Spin for a few seconds in a microcentrifuge to bring the contents to the bottom of the tube.

8. Incubate at 37°C for 10 minutes.

7. Avoid vigorous mixing of the reaction mixture as this can cause severe loss of enzyme activity.

8. Purified DNA can be labelled to high specific activity in 10 minutes at 37°C but, if desired can be labelled for up to 1 hour at this temperature.

Protocol

Notes

8. Incubate at 37°C for 10 minutes *continued*.
9. Stop the reaction by the addition of 5 µl of 0.2 M EDTA. For use in a hybridization, denature the labelled DNA by heating to 95–100°C for 5 minutes, then chill on ice.

8. **Continued**
When labelling DNA in low melting point agarose, longer incubation of 15–30 minutes at 37°C are required for optimum labelling. Longer incubation times (up to 60 minutes) are required when nucleotide analogues (e.g. [³⁵S]dNTP(S)) are used.
9. Labelled probe can be stored at -15°C to -30°C in a non frost-free freezer. Prolonged storage of ³²P-labelled probes can lead to substantial probe degradation(7). High specific activity probes should be stored for no longer than 3 days. Although probe purification is not usually necessary for most membrane applications the removal of unincorporated nucleotide is sometimes useful to reduce background in filter hybridizations for probes >10⁹ dpm/µg or when the reaction yields an incorporation of less than 50%. This procedure is

Protocol

Notes

9. Stop the reaction by the addition of 5 μ l of 0.2 M EDTA. For use in a hybridization, denature the labelled DNA by heating to 95-100°C for 5 minutes, then chill on ice *continued*.

9. *Continued*

described in Appendix III. Calculation of probe specific activity is described in Appendix II. Extensive experimentation with **Rapid-hyb buffer (RPN1635/6)** has shown that probe purification, even under the conditions given above is not required with the isotopes ^{32}P and ^{32}P . Purification of ^{32}S labelled probes is however required to reduce filter background.




Figure 2. Time course of incorporation of [α - 32 P]dCTP (17 pmol) in a Megaprime reaction at 37°C. The DNA used was the standard DNA supplied with the system.

5.3. Use of alternative reaction conditions

a. Use of more than one labelled [α - 32 P]dNTP.

Table 1 lists the results of a selection of standard reactions, using a variety of input labels under optimum conditions. Figure 3 gives more complete information on their use in Megaprime reactions. Reactions were carried out at 37°C for 5 minutes.

b. Use of alternative specific activity [α - 32 P]dNTPs.

When using [α - 32 P]dNTPs of specific activity <3000 Ci/mmol the incubation time should be extended to 1 hour at 37°C.

c. Use of [³²P]dNTP α S.

When using ³²S-labelled radionucleotides the incubation time should be extended to 1 hour at 37°C.

d. Labelling at room temperature.

If desired, labelling reactions can be carried out at room temperature. Maximum incorporation occurs after an incubation time of 45–60 minutes. A decline in incorporation can be observed if reactions are left overnight.

e. Factors affecting the labelled DNA.

1. Specific activity

Figure 3a should be used to ascertain the number and quantity of labelled dNTP's required in order to prepare a probe of the desired specific activity.

2. Efficiency

Figure 3b indicates the efficiency of the chosen reaction conditions, and thus permits a balance of specific activity and economy.

3. Probe length

Figure 3c gives a measure of mean probe lengths obtained under standard conditions. Probe lengths were measured by denaturing agarose gel electrophoresis followed by autoradiography with reference to molecular weight standards.

Probe length can be affected by the concentration of DNA, primer and nucleotide, the size of the template DNA and also radiolysis of the labelled probe. The data in the figure was obtained using linearized plasmid DNA, 4.5 Kb in length under the standard labelling conditions.

It is recommended that not less than 10 pmol and not more than 125 pmol of any labelled dNTP is used in the reaction and combinations shown offer optimum balance of stability, specific activity and economy.

Compounds	Specific activity (see note a)	Formulation (see note b)	Product code	Quantity of each dNTP required	Specific activity of probe (see notes c,d and e)	
	TBq/mmol			MBq		
	Ci/mmol			μCi		
				μl		
				pmol		
[α- ³² P]dCTP	~111 ~222	~3000 ~6000	1 1	PB 10205 PB 10475	1.85 50 5 17 7.4 200 20 32	1.9 × 10 ⁹ 5.3 × 10 ⁹
[α- ³² P]dATP	~111	~3000	1	PB 10204	1.85 50 5 17	3.4 × 10 ⁹
[α- ³² P]dCTP	~111	~3000	1	PB 10205	1.85 50 5 17	
[α- ³² P]dATP	~111	~3000	1	PB 10204	1.85 50 5 17	3.7 × 10 ⁹
[α- ³² P]dCTP	~111	~3000	1	PB 10205	1.85 50 5 17	
[α- ³² P]dGTP	~111	~3000	1	PB 10206	1.85 50 5 17	

a. At the specific activity reference date of the labelled nucleotide.


b. Formulation code 1 = 370 MBq/ml, 10 mCi/ml in stabilized aqueous solution.

c. The probe specific activities were calculated using observed incorporation levels which are similar to those found in figure 3b.

d. It is important to note that the specific activity of probes made from different amounts of labelled dNTP cannot be calculated on a proportional basis, because net DNA synthesis occurs.

e. Brackets enclose nucleotides used in combination.


a) Specific activity



- i) One labelled dNTP
- ii) Two labelled dNTP
- iii) Three labelled dNTP

Figure 3. The use of [α -³²P]dNTPs in the Megaprime DNA labelling system (see notes on page 26).


b) Incorporation efficiency



- i) One labelled dNTP
- ii) Two labelled dNTP
- iii) Three labelled dNTP

Figure 3. The use of $[\alpha\text{-}^{32}\text{P}]\text{dNTPs}$ in the Megaprime DNA labelling system (see notes on page 26).

c) Probe length



- i) One labelled dNTP
- ii) Two labelled dNTP
- iii) Three labelled dNTP

Figure 3. The use of [α - 32 P]dNTPs in the Megaprime DNA labelling system (see below).

Notes to figure 3

- a.** The results shown are the means of a number of experiments in which different nucleotides and combinations of nucleotides were used. Observed results may deviate $\pm 10\%$ from those shown.
- b)** As the number of different labelled nucleotides is increased, at a given level of total input label, the net synthesis of DNA is reduced. Although the overall incorporation efficiency is reduced the labelled product is of a higher specific activity.

c. The data was generated using the standard labelling protocols.

If dNTPs <3000 Ci/mmol are to be used, then the desired probe specific activity must be multiplied by a conversion factor, before determining the amount of input label.

For a single labelled dNTP:-

$$\text{Total input label (pmols)} = \frac{3000 \text{ Ci/mmol}}{\text{specific activity of dNTP to be used}} \times \text{required probe specific activity}$$

For more than one labelled dNTP the mean specific activity of the labelled dNTP to be used should be inserted in the above calculation.

Having determined the required number of pmols of input label with reference to figure 3a, the required volume of each labelled dNTP can be calculated. Note that the figures give the **total** amount of input label required. If more than one labelled dNTP is

to be used, this figure should be divided by the number of labelled dNTPs to be used to give the required number of pmols of each labelled dNTP.

Volume of each labelled dNTP required in μl =

$$\frac{\text{pmol of dNTP required} \times \text{specific activity of dNTP (Ci/mmol)} \times 10^{-3}}{\text{radioactive concentration of dNTP (mCi/ml)}}$$

6. Appendices

6.1. Appendix I. Labelling of DNA fragments in low melting point agarose

The DNA samples produced by the following protocol have been found to be labelled to approximately the same extent as purified DNA. 15–20 minutes at 37°C is optimum for labelling. The standard labelling protocol may be found to be more appropriate for labelling DNA in agarose as the volume of DNA to be added using the new protocol is limited to 5 µl, requiring a relatively high initial DNA concentration.

Protocol	Notes
<p>1. Fractionate restriction endonuclease digested DNA in a suitable low melting point agarose gel containing 0.5 µg/ml ethidium bromide. Estimate the DNA content of the band by reference to a set of standards of known concentration on another track. 250 ng should allow 25 ng to be used in the standard labelling protocol without further concentration</p>	<p>1. A low melting point agarose of high purity for example SepRate-LMP is recommended for maximum labelling efficiency.</p>
<p>2. Excise the desired band cleanly, with the minimum of excess agarose and transfer to a pre-weighed 1.5 ml microcentrifuge tube.</p>	<p>2. It is recommended that the exposure to UV light is minimized, as prolonged exposure can damage the DNA.</p>

Protocol	Notes
<p>3. Add water to a ratio of 3 ml per gram of gel and place in a boiling water bath for 5 minutes to melt the gel and denature the DNA.</p>	<p>3. If the DNA is not to be used immediately divide the boiled samples into suitably sized aliquots and store at -15°C to -30°C in a non frost-free freezer.</p>
<p>4. If the DNA is to be used immediately remove the appropriate volume containing 25 ng, add to the primers as indicated in the labelling protocol (page 11, step 3). The volume of DNA should not exceed 25 µl for the standard labelling protocol.</p>	<p>4. When using DNA which has been previously boiled and then stored at -15°C to -30°C, first place the tube in a boiling water bath for 30 seconds to melt the agarose, before removing the required volume containing 25 ng. Do not reboil DNA aliquots more than twice.</p>
<p>5. Incubate the labelling reaction for 15–20 minutes at 37°C.</p>	

6.2. Appendix II. Monitoring the reaction and calculating the specific activity of the labelled DNA

A. Adsorption to DE81 paper

Monitoring of the progress of the labelling reaction and measurement of probe specific activity can be achieved by determining the proportion of the radionucleotide incorporated during the Megaprime reaction.

1. Remove a 1 or 2 μl aliquot of the reaction mixture to a clean microcentrifuge tube containing 20 μl of water or 10 mM Tris/HCl pH.8.0. 1 mM EDTA buffer. Mix well by pipetting up and down.
 2. Spot, in quadruplicate, 5 μl aliquots of this dilution on to Whatman DE81 chromatography paper squares (minimum size 1 x 1 cm), placed on a non-absorbent backing. These squares may be marked with a pencil for identification if required.
 3. Take two of the filters and dry under a heat lamp. 10–15 minutes should be adequate.
 4. Wash the remaining two filters twice for 5 minutes each, at room temperature in excess 2xSSC (30 mM $\text{Na}_3\text{citrate}$, 300 mM NaCl pH7.0) using gentle agitation. Rinse briefly in distilled water and then once with ethanol for 5 minutes. Then dry the filters under a heat lamp.
4. In aqueous solution DE81 paper becomes fragile and care should be taken when handling. In order to stabilize the paper the squares are rinsed in ethanol.

Protocol

5. Place the squares in separate vials with at least 5 ml of scintillation fluid and count.
6. Efficiency of counting will vary, but the percentage incorporation can be used to calculate probe specific activity. Unlike the nick translation labelling reaction, Megaprime labelling leads to net DNA synthesis, and so the total amount of DNA at the end of the reaction must be calculated.

Total amount of DNA (A) ng =

$$\frac{\text{Total number of } \mu\text{Ci added} \times 13.2 \times \% \text{ incorporation} + 25}{\text{Number of radioactive dNTPs added} \times \text{average specific activity of dNTPs added}}$$

Number of radioactive dNTPs added x average specific activity of dNTPs added

This assumes a 25% content of any one dNTP in the newly synthesized DNA, and 25 ng of template DNA.

*13.2 equals four times the average molecular weight of the four dNTPs divided by 100.

Notes

5. Determination of the proportion of the ^{32}P labelled nucleotide incorporated may be achieved using Cerenkov counting if desired in this case drying the filter is not necessary.
6. The mean value of the counts on the washed filter represents the proportion of the radionucleotide incorporated into the DNA probe, while the mean of the unwashed filters represents the total amount of radioactivity in the reaction mix, such that;

$$\% \text{ incorporation} = \frac{\text{mean counts on washed filters} \times 100}{\text{mean counts on unwashed filters}}$$

6. Continued.

The amount of radioactivity incorporated during the reaction (B) in dpm.

$B = \text{total number of } \mu\text{Ci added} \times 2.2 \times 10^4 \times \% \text{ incorporation}$

Thus the specific activity of the labelled DNA is

specific activity = $\frac{B}{A} \times 10^3$ dpm per μg

B. Precipitation with trichloroacetic acid

Plastic or siliconized glass tubes must be used to avoid adsorption of DNA.

1. Dilute an appropriate aliquot of the reaction mixture as described in section A1.
2. Transfer 1–10 μl of diluted reaction mixture to two duplicate tubes containing 200 μl water or 0.2M EDTA and 50 μl carrier DNA solution. Mix well. Use this mixture (less any set aside in step 3) for the TCA precipitation described in step 4 below.
3. Set aside an appropriate aliquot from each tube in step 2 for the determination of total input radioactivity.
4. To the diluted samples from step 2, add 2 ml ice-cold 10% trichloroacetic acid (TCA) solution, vortex, and allow to stand in an ice-bath for 10–15 minutes. The labelled and carrier DNA will co-precipitate. **Note that TCA is corrosive, and care should be taken in its handling.**
5. Collect the precipitated DNA by vacuum filtration on a glass fibre or nitrocellulose filter disc.

6. Wash the filter discs six times with 2 ml 10% TCA solution and dry the filter discs thoroughly, for example using an infra-red lamp. Avoid overheating and possible charring of the discs.
7. Count the dried filter discs by liquid scintillation or Cerenkov (^{32}P) and count with the samples set aside in step 3.
8. Determine % incorporation and probe specific activity as in section A6.

6.3. Appendix III. Removal of unincorporated nucleotides

Removal of unincorporated nucleotides is sometimes desirable to reduce background produced by the probe during hybridization. It is considered important to remove these free nucleotides particularly if the radioactive probe is to be kept for several days before use or the incorporation is less than 50%. If ^{32}P or ^{33}P -labelled probes are to be used in combination with GE Healthcare's new Rapid-hyb buffer (RPN1635/6), purification is not required unless the probe is to be used more than 24 hours after preparation. Probes can be purified by Sephadex chromatography or selective precipitation (8,9).

A. Sephadex™G-50 spin columns

Probe reaction are passed through columns packed with Sephadex G-50, which retains the free nucleotides within the column matrix. A number of pre-packed columns are commercially available. However columns may also be prepared as indicated below:

1. Equilibrate Sephadex G-50 in TE buffer either overnight or at 65°C for 1–2 hours.
2. Plug a 1.0 ml syringe with a piece of siliconized glass wool.
3. Fill the syringe with the equilibrated Sephadex. Place in a 15 ml conical tube, in which a decapped 1.5 ml microcentrifuge tube has been inserted. Centrifuge at 1600 g for 5 minutes. Remove

any liquid from the microcentrifuge tube. Refill with Sephadex and centrifuge as before. Continue until the column is packed to a volume of 1 ml.

TM Sephadex is a trademark of GE Healthcare

4. Add a volume of TE buffer equal to the reaction volume, to the top of the column and centrifuge, as in step 3. A minimum of 50 μ l should be applied to the column.
5. Repeat once more to ensure fractions of the correct size are collected from the column.
6. Place the column in a clean 15 ml conical tube containing a decapped 1.5 ml microcentrifuge tube.
7. Apply the DNA sample to the column. Centrifuge as before. The purified probe is collected in the microcentrifuge tube.

B. Selective precipitation of labelled DNA

The following protocol leads to precipitation of DNA greater than about 20 nucleotides in length with unincorporated nucleotides remaining in solution. Recovery of the labelled DNA by this method varies according to the DNA concentration and size, and may be as low as 50%.

1. Add one volume of 4 M ammonium acetate, pH4.5 to the nick translation reaction, and mix gently by pipetting up and down.
2. Add four volumes of ethanol, mix by inversion. Chill the mixture for 15 minutes in a dry-ice ethanol bath or place at -70°C for at least 30 minutes.
3. Thaw the mixture if necessary by placing at 37°C for 2 minutes.
4. Spin in a microcentrifuge for 15 minutes. Carefully aspirate and dispose of supernatant in a suitable manner.
5. Wash the pellet once in 0.5 ml of 0.67 M ammonium acetate, pH 4.5, 67% ethanol at room temperature by gentle inversion, centrifugation and aspiration.

6. Wash the pellet once in 90% ethanol, in the same manner. Dry the pellet.
7. Finally redissolve the DNA pellet in TE buffer for use as a probe and for storage.

6.4. Appendix IV. Additional equipment and reagents

TE buffer (10 mM Tris/HCl, pH 8.0, 1 mM EDTA)

0.2 M EDTA solution

Adjustable pipettes for example Pipetman™

Sterile pipette tips

Waterbaths at 37°C and 100°C

Polypropylene microcentrifuge tubes

Microcentrifuge

Gloves

Radiation safety equipment

DE81 ion-exchange chromatography paper (Whatman)

Trichloroacetic acid (TCA) solution: 10% (w/v) TCA in water

Filter discs; glass fibre or nitrocellulose

Plastic or siliconized glass tubes, capacity ~5 ml

Filtration apparatus

2x SSC (30 mM Na₃ citrate, 300 mM NaCl, pH 7.0)

™ Pipetman is a registered trademark of Gilson

7. Troubleshooting guide

If poor results are obtained, the following guide may help to determine the cause of the problem.

Problem	Possible cause	Remedy
1. Low signal	1. Incomplete denaturation of template DNA	1. Ensure denaturation protocol is followed.
	2. Low probe concentration	2. Accurately measure the concentration of template DNA used in the labelling reactions. Check recovery of probe if purification is performed to remove unincorporated nucleotide.
	3. Low probe specific activity	3. If the specific activity of the labelled DNA is lower than expected, a labelling reaction should be carried out using a sample of the control DNA supplied with the system. If this proceeds satisfactorily, check the concentration and purity of your DNA.

Problem	Possible cause	Remedy
	4. Loss of dNTP during evaporation	4. If the dNTP solution has been evaporated to dryness prior to use, handling losses may have occurred. Check this loss has not occurred during lyophilization of the solvent, during transfer of the reconstituted dNTP solution or by adsorption of the dNTP onto the walls of the tube. If necessary the reconstituted dNTP solution may be counted and an adjustment made before setting up the labelling reaction.
2. Non-specific background over whole of filter	1. Presence of unincorporated label	1. Unincorporated nucleotides can give high backgrounds. Remove by Sephadex G-50 spin columns or ethanol precipitation (see page 32 for protocol)

Problem	Possible cause	Remedy
	2. Concentrated probe has contacted membrane directly during probe addition	2. It is suggested that up to 1.0 ml of the buffer used for prehybridization is withdrawn for mixing with the probe. The mixture should then be added back to the hybridization container in an area away from the filter.
	3. Probe concentration is too high	3. Ensure measurement of template DNA concentration is accurate
	4. Probe not denatured	4. Non-denatured double-stranded probes often give high backgrounds.

8. References

1. FEINBERG, A.P. and VOGELSTEIN, B., *Anal. Biochem.*, **132**, pp.6-13, 1983.
2. FEINBERG, A.P. and VOGELSTEIN, B., *Addendum Anal. Biochem.*, **137**, pp.266-267, 1984.
3. SOUTHERN, E.M., *J.Mol.Biol.*, **98**, pp.503-517, 1975.
4. THOMAS, P.S., *Proc. Natl. Acad. Sci. USA.*, **77**, pp.5201-5205, 1980.
5. MEINKOTH, J. and WAHL, G., *Anal. Biochem.*, **138**, pp. 267-284, 1984.
6. GRUNSTEIN M. and HOGNESS, D.S., *Proc. Natl. Acad. Sci. USA.*, **72**, pp. 3961-3965, 1975.
7. HODGSON, C.P., FISK, R.Z. and WILLET, L.B., *Biotechniques*, **6**, pp.208-211.
8. SAMBROOK, J. FRITSCH, E.F. and MANIATIS, T., *Molecular Cloning*, a laboratory manual (second edition), Cold Spring Harbour Laboratory, 1989.
9. MUNDY, C.R., CUNNINGHAM, M.W. and READ, C.A., *Essential Molecular Biology; A Practical Approach Vol 2* (T.A. Brown, ed) Oxford University Press, Oxford, 1991. pp.57-109.

9. Related Products

Labelling systems

Nick translation kits	N5000/5500
3'-end labelling kit	N4020
5' end labelling kit	RPN 1509
RNA labelling system (paired promoter SP6/T7 system)	RPN 3100

Hybridization buffers

Rapid-hyb buffer	RPN 1635/6
Hybridization buffer tablets	RPN 131

Hybridization membranes

Hybond™ - Range of nylon and nitrocellulose blotting membranes

Autoradiography products

Hyperfilm™ - high performance autoradiography films
Hypercassettes™ and Hyperscreens™ - available from stock

Safety Products

Radiation safety products for safe handling and storage of $^{32}\text{P}/^{33}\text{P}/^{35}\text{S}$ and ^{125}I , liquid scintillation products

Agarose

SepRate™ - range of highly purified agarose for a range of DNA fragment sizes and users

Labelled dNTPs

See Table 2

Table 2. Labelled dNTPs and analogues available from GE Healthcare

Compound	Specific Activity		Formulation (see key)	Product code
	TBq/mmol	Ci/mmol		
[α - ³² P]dATP	~220	~6000	1	PB 1074
	~110	~3000	1	PB 10204
			2	PB 204
	~30	~800	1	PB 10384
	~15	~400	1	PB 10164
			2	PB 164
[α - ³² P]dCTP	~220	~6000	1	PB 10475
	~110	~3000	2	PB 10205
			1	PB 205
	~30	~6000	1	PB 10385
	~15	~400	2	PB 10165
				PB 165
[α - ³² P]dGTP	~220	~3000	1	PB 10206
			2	PB 206
			1	PB 10386
	~30	~800	1	PB 10386
	~15	~400	1	PB 10166
			2	PB 166
[α - ³² P]dTTP	~110	~300	1	PB 1027
			2	PB 207
			1	PB 10387
	~30	~800	1	PB 10387
	~15	~400	1	PB 10167
			2	PB 167
[α - ³² P] dATP	37-110	1000-3000	1	BF 1001
[³⁵ S]dATP α S	>37	>1000	1	SJ 1304
	~22	~600	1	SJ 304
	~15	~400	1	SJ 264

Compound	Specific Activity		Formulation (see key)	Product code
	TBq/mmol	Ci/mmol		
[³⁵ S]dCTPαS	>37	>1000	1	SJ1305
	~22	~600	1	SJ 305
	~15	~400	1	SJ 265
[³⁵ S]dGTPαS	~22	~600	1	SJ 306
[³⁵ S]dTTPαS	~22	~600	1	SJ 307
[8- ³ H]dATP	0.37-1.1	10-30	2	TRK 347
[1',2',2,8- ³ H]dATP	1.83-3.7	50-100	2	TRK 633
[1',2',5- ³ H]dCTP	1.85-3.14	50-85	2	TRK 625
[5- ³ H]dCTP	0.55-1.1	15-30	2	TRK 352
[8- ³ H]dGTP	0.185-0.740	5-20	2	TRK 350
[1'2,(- ³ H]dGTP	0.9-1.85	25-50	2	TRK 627
[methyl, 1',2'- ³ H]TTP	3.3-4.8	90-130	2	TRK 576
[methyl- ³ H]TTP		40-60	2	TRK 424
		30	2	TRK 354
[¹²⁵ I]dCTP	>55	>1500	3	IM 5103

Formulation codes:

1) = 370 MBq/ml, 10 mCi/ml in stabilized solution

2) = 37 MBq/ml, 1 mCi/ml in 50% aqueous ethanol

3) = 185 MBq/ml, 5 mCi/ml in 50% aqueous ethanol

See GE Healthcare Products catalogue for further details.

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30 751 84

Uppsala

Sweden

GE Healthcare Europe GmbH

Munzinger Strasse 5 D-79111

Freiburg

Germany

GE Healthcare UK Limited

Amersham Place

Little Chalfont

Buckinghamshire

HP7 9NA

UK

GE Healthcare Bio-Sciences

Corp

800 Centennial Avenue

P.O. Box 1327

Piscataway

NJ 08855-1327

USA

GE Healthcare Bio-Sciences KK

Sanken Bldg. 3-25-1

Hyakunincho Shinjuku-ku

Tokyo 169-0073

Japan

GE Healthcare regional office contact numbers:**Asia Pacific**

Tel: +85 65 62751830

Fax: +85 65 62751829

Australasia

Tel: + 61 2 8820 8299

Fax: +61 2 8820 8200

Austria

Tel: 01/57606-1613

Fax: 01/57606-1614

Belgium

Tel: 0800 73 890

Fax: 02 416 8206

Canada

Tel: 1 800 463 5800

Fax: 1 800 567 1008

Central, East, & South**East Europe**

Tel: +43 1 972 720

Fax: +43 1 972 722 750

Denmark

Tel: 45 70 25 24 50

Fax: 45 45 16 2424

Eire

Tel: 1 800 709992

Fax: +44 1494 542010

Finland & Baltics

Tel: +358 9 512 3940

Fax: +358 9 512 39439

France

Tel: 01 69 35 67 00

Fax: 01 69 41 98 77

Germany

Tel: 0800 9080 711

Fax: 0800 9080 712

Greater China

Tel: +852 2100 6300

Fax: +852 2100 6338

Italy

Tel: 02 26001 320

Fax: 02 26001 399

Japan

Tel: +81 3 5331 9336

Fax: +81 3 5331 9370

Korea

Tel: 82 2 6201 3700

Fax: 82 2 6201 3803

Latin America

Tel: +55 11 3933 7300

Fax: + 55 11 3933 7304

Middle East & Africa

Tel: +30 210 96 00 687

Fax: +30 210 96 00 693

Netherlands

Tel: 0800-82 82 82 1

Fax: 0800-82 82 82 4

Norway

Tel: +47 815 65 777

Fax: +47 815 65 666

Portugal

Tel: 21 417 7035

Fax: 21 417 3184

Russia, C.I.S. & N.I.S

Tel: +7 495 956 5177

Fax: +7 495 956 5176

Spain

Tel: 902 11 72 65

Fax: 935 94 49 65

Sweden

Tel: 018 612 1900

Fax: 018 612 1910

Switzerland

Tel: 0848 8028 10

Fax: 0848 8028 11

UK

Tel: 0800 515 313

Fax: 0800 616 927

USA

Tel: +1 800 526 3593

Fax: +1 877 295 8102

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited

Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA

UK



imagination at work