

Sequenase Version 2.0 DNA Sequencing Kit

Product Number 70770

Brief Protocol

1. Denature double-stranded templates.

2. **Annealing reaction:**

DNA (approx. 0.5pmol)	__μl (up to 7μl)
H ₂ O	__μl (to adjust total volume)
Sequenase™ Reaction Buffer	2μl
Primer (0.5pmol)	1μl
Total	10μl

Anneal by heating 2 minutes at 65°C, then cool slowly to <35°C over 15-30 minutes. Chill on ice for use in step 6.

3. While cooling, label, fill and cap tubes with 2.5μl of each termination mix. Use mixes from **red-capped** tubes for dGTP or **orange-capped** tubes for dTTP. Keep covered at room temperature for steps 5 and 7.

G (2.5μl)	A (2.5μl)
T (2.5μl)	C (2.5μl)

4. Dilute labeling mix 1:5 to working concentration, dGTP (**green-capped** tube) or dTTP (**yellow-capped** tube).

Retain for use in step 6.

Labeling mix	__μl (typically 2μl)
H ₂ O	__μl (typically 8μl)

5. Pre-warm 4 termination tubes from step 3 ('G', 'A', 'T' and 'C') in a 37°C bath.

6. **Labeling reaction**

To ice-cold annealed DNA mixture (10μl), add:

Dithiothreitol (DTT), 0.1M	1μl
Diluted labeling mix	2μl
[α- ³⁵ S or α- ³³ P]dATP	0.5μl
Diluted Sequenase polymerase*	2μl
Total	15.5μl

Mix and incubate at room temperature 2-5 minutes.

7. **Termination reactions**

Transfer 3.5μl of the labeling reaction to each termination tube ('G', 'A', 'T' and 'C'), mix and continue incubation of the termination reactions at 37°C for 5 minutes.

8. Stop the reactions by adding 4μl of stop solution.

9. Heat samples to 75°C for 2 minutes immediately before loading onto sequencing gel.

* **Important:** The Sequenase DNA polymerase must be diluted 8-fold prior to use. If desired, the polymerase for the complete kit can be pre-diluted as follows:

Sequenase DNA polymerase (70775)	25μl
Pyrophosphatase (70950)	25μl
Glycerol Enzyme Dilution Buffer (70799)	150μl

Using this dilution will necessitate the use of a Glycerol Tolerant Sequencing Gel (see protocol book). The polymerase may also be diluted as needed by mixing 1μl polymerase with 0.5μl pyrophosphatase and 6.5μl of enzyme dilution buffer (70765). Polymerase diluted this way may be used for ordinary Tris-Borate-EDTA (TBE) gels.

All goods and services are sold subject to the terms and conditions of sale by USB Corporation and the group which supplies them. A copy of these terms and conditions is available on request.

Refer to Protocol Booklet, 70770A, for detailed information on licensing and patents.

Sequenase is a trademark of GE Healthcare Bio-Sciences Ltd or its subsidiaries.

USB and the logo design are registered trademarks of USB Corporation.

©USB Corporation 2000 - All rights reserved

Printed in the United States

USB Corporation
26111 Miles Road
Cleveland, Ohio 44128 USA
www.usbweb.com



GE imagination at work



96-0770B/00/10